



HEARING MATTERS AUSTRALIA

Support from real life experience

Ground Floor, Suite 600
The Australian Hearing Hub
16 University Ave
Macquarie University NSW 2109
(02) 9878 089
admin@hearingmattersaustralia.org
www.hearingmattersaustralia.org

BLUETOOTH COMPATIBLE HEARING AIDS

What is Bluetooth?

Bluetooth technology was developed by Ericsson, for use with mobile phones. Bluetooth uses radio frequencies to allow devices such as mobile phones, digital tablets, computers and smart television sets to communicate with each other without the use of any cable or wiring. Typically, Bluetooth technology can be used over a maximum range of 10 metres. Bluetooth creates a local wireless network between multiple electronic devices for multiple purposes. Initially the two devices that want to communicate are paired, thus creating a secure radio link that is encrypted to prevent other devices from interfering.



How does a Bluetooth hearing aid work?

Bluetooth hearing aids have Bluetooth wireless technology built into them and it cannot be added afterwards. A small communication device, usually worn around the neck on a neck loop, is also sometimes required as an interface between the hearing aid and the other device. There are two ways to connect to other devices: one-way communication such as listening to the radio or television and two-way communication such as using a telephone. Two-way communication uses a microphone that is built into the communication device on the neck loop or into the hearing aid itself (depending on the model).

How do you use a Bluetooth hearing aid?

The procedure for 'pairing' the other Bluetooth compatible device with the hearing aid will vary depending on the manufacturer. Once the devices are 'paired' they should not need pairing again. When connected to a Bluetooth compatible phone, it will ring inside the hearing aids providing you are in range. If two hearing aids are worn, you will hear the phone call in both ears. When connected to a one-way audio device, an additional Bluetooth transmitter is required to connect the other device such as a television or iPod to the hearing aid via Bluetooth.

Information sheets are shared by HMA as a service to members and those in the community who have an interest in hearing loss. Every effort has been made to ensure the accuracy of the information provided, however HMA accepts no responsibility for any adverse consequences arising from the contents of these sheets.

HMA information sheets are for personal use only. Downloads and printing allowable for whole sheets.

Contact HMA with queries or amendments.

Find out more and join HMA for additional member benefits at www.hearingmattersaustralia.org HMA INFO SHEET 10/2020



What can Bluetooth be used for?

Devices that can be used with a Bluetooth compatible hearing aid, providing they are also Bluetooth compatible, include:

- Mobile Phone
- Bluetooth Home Phone
- Television
- Radio
- Computer
- -iPod or Android tablet
- -MP3 Player

A Bluetooth hearing aid probably also has a T-Switch, which enables it to be used with a loop or in locations that are equipped for hearing aid use.

Who is a Bluetooth hearing aid suitable for?

Bluetooth compatible hearing aids are not suitable for everyone. An audiologist will be able to describe the benefits of Bluetooth hearing aids for the individual and may be able to provide a demonstration of how the technology is used. Not everyone will find the Bluetooth technology valuable, but many enjoy the ease of use of Bluetooth technology for their hearing aids.

What should I think about before I buy a hearing aid with Bluetooth capabilities?

Before purchasing a Bluetooth compatible hearing aid, talk to your audiologist to discuss whether it is a good option for you. Consider your lifestyle and what devices you would use the technology with. Keep in mind that as well as the hearing aids, there are other accessories to keep track of, and for which you will need to charge batteries. You will also need to check that other gadgets such as television and radio have Bluetooth capabilities, as some older models may not.